

EasyTEMP Free vs. Professional

	EasyTEMP Free	EasyTEMP Professional
Control of one JULABO unit with integrated interface	✓	✓
Control of up to 24 JULABO units with integrated interface		✓
Instrument:		
Recording and display of currently measured values on PC	✓	✓
Setpoint programming via PC	✓	✓
Status Indication	✓	✓
Individual control window for each unit		✓
Simultaneous start of units with just one button		✓
Recording of Measured Values:		
Graphic zoom data function	✓	✓
Displays up to 4 curves in one diagram	✓	
Displays unlimited curves in one diagram		✓
Curves can be assigned to individual scales		✓
Insert formulas such as averaging, differences between measured values, etc.		✓
Insert text comments with display in the diagram		✓
Ramp Programming:		
Ramp function up to 100 steps	✓	
Ramp function up to 1000 steps (individually for each instrument)		✓
Cyclic repetition of stored profiles	✓	✓
Modify running profiles	✓	✓
Graphic display of total profiles		✓
Data Recording:		
Records measured values in ASCII format	✓	✓
Records measured values in Microsoft® Excel®		✓
Saves additional, relevant measuring data		✓
Scalable recording pattern		✓
Export function of graphs into JPG format		✓
Uploads previously created recordings with print function		✓
Control and integration of laboratory instruments of different makes e.g. stirrers, balances, dosing pumps, pH meters, etc.		✓