

CORIO C-BT27


The open heating bath circulators of the new CORIO series feature highquality transparent bath tanks.

Product features

- Models for internal applications
- Bright, white, easy to read display
- Very quiet
- Precise temperature control



Heat-up time (230V, 2kW heater)



Performance values

115V/60Hz (Nema N5-15 Plug)	
Heating capacity kW	1
Viscosity max. cST	1
Circulation capacity flow rate l/min	6
Circulation capacity pressure psi	1.5
Power A	10

Technical data

Available voltage versions		Bath	
Order No.	9 011 327	Bath tank	Polycarbonat
Available voltage versions:		Usable bath opening in. (W x L / D)	11.8 x 13.8 / 7.9
9011327.22	100-115V/50-60Hz (Nema N5-15 Plug)		



9011327.33.chn	230V/50-60Hz (CN Plug)
9011327.33	230V/50-60Hz (Schuko Plug - CEE 7/4 Plug Type F)
9011327.04	230V/50-60Hz (UK Plug Type BS1363A)
9011327.05	230V/50-60Hz (CH Plug Type SEV 1011)

Other		Electronics	
Classification	Classification I (NFL)	Temperature control	PID1
IP Code	IP 21	Absolute temperature calibration	1 Point Calibration
Pump type	Immersion Pump	Temperature displayTemperature display	LED
		Temperature settingTemperature setting	Keypad
Dimensions and volumes		Temperature values	
Weight lbs	15.9	Setting the resolution of the temperature display $^{\circ}\mathrm{C}$	0.1
Dimensions in. ($W \times L \times H$)	15 x 22.8 x 16.9		
Filling volume I 20	20 27	Working temperature range °C	+20 +100
	20 27	Temperature stability °C	±0.03
		Ambient temperature °C	+5 +40
		Temperature display resolution °C	0.01 0.1

For applications near or below ambient temperature: use a cooling coil or JULABO immersion cooler.

All Benefits



100% Checked.

100% testing. 100% quality. Each JULABO Circulator undergoes thorough quality testing before leaving the factory.



JULABO. Quality. Highest standards of quality for a long product life.



Satisfied customers.

11 subsidiaries and more than 100 partners worldwide guarantee fast and qualified JULABO support.



Stable.

Rubber feet allow for a secured footing of your CORIO to prevent damage to your laboratory equipment.



Brilliant. Very bright display makes it easy to read even from a distance.





Development consistently applied environmentally friendly materials and technologies.



Modern. Reliable.

High-grade components in every JULABO Circulator - platinum sensors, proven motor technology, CAN-Bus communication and much more.



Services 24/7.

Around the clock availability. You can find suitable accessories, data sheets, manuals, case studies, and more at www.julabo.us.

Quick start.

Individual JULABO consultation and comprehensive manuals at your disposal.

Connection. Easy.



Inclined pump connections (M16×1) facilitate the connection of applications. Each unit includes 2 barbed fittings of 8/12 mm diameter each.





Everything at the front.

All operating controls and safety functions are accessed easily and comfortably from the front.



Locked in. The lockable power plug guarantees safe

connection. More process safety.



Timer. Integrated.

CORIO circulators include an integrated timer function. When the set time has elapsed, a signal sounds and the device switches off. Setting range: 0 ... 999 minutes.



Exact.

You can rely on it. PID1 control and 'Active Cooling Control' make the new CORIO precise and perfect.



Switch on. And off you go. Intelligent operating concept. Ready for operation with just a few quick and easy steps.



Early warning system for low liquid level Maximum safety for applications, optical and audible alarm, allows user to refill bath fluid before the unit shuts down