Julabo

CORIO CD Heating immersion circulator

In addition to the advantages of the CORIO C circulator, the CORIO CD immersion circulator can be equipped with an optional pump set for temperature control of external applications.

Product features

- Precise temperature control
- For internal and external applications (accessories required)
- For bath tanks with a max. filling volume of 50 liters
- Bright, white, easy to read display
- Very quiet
- Easy operation
- Easy pump change-over between internal and external circulation
- USB connection
- Class III (FL) according to DIN 12876-1
- Powerful and infinitely adjustable pressure pump



Performance values

115V/60Hz (Nema N5-15 Plug)	
Heating capacity kW	1
Viscosity max. cST	50
Pump capacity flow rate l/min	16
Pump capacity flow pressure psi	4.8
Power A	10

230V/60Hz (Schuko Plug - CEE 7/4 Plug Type F)	
Heating capacity kW	2
Viscosity max. cST	50
Pump capacity flow rate I/min	17
Pump capacity flow pressure psi	6.2
Power A	10

Technical data

Available voltage ve	rsions	Cooling	
Order No.	9 012 000	Cooling coil	accessory
Available voltage version	s:		
9012000.01	100V/50-60Hz (Nema N5-15 Plug)		
9012000.02	115V/60Hz (Nema N5-15 Plug)		
9012000.03	230V/50Hz (Schuko Plug - CEE 7/4 Plug Type F)		
9012000.03.chn	230V/50Hz (CN Plug)		
9012000.04	230V/50Hz (UK Plug Type BS1363A)		



9012000.05	230V/50Hz (CH Plug Type SEV 1011)
9012000.13	230V/60Hz (Schuko Plug - CEE 7/4 Plug Type F)

Other		Electronics	
Classification	Classification III (FL)	Digital interface	USB
IP Code	IP 21	Temperature control	PID1
Pump function	Pressure Pump	Absolute temperature calibration	1 Point Calibration
Pump type	Immersion Pump	Temperature displayTemperature display	LED
		Temperature settingTemperature setting	Keypad
Dimensions and volumes		Temperature values	
Dimensions and volumes Usable immersion depth in.	6.5	Setting the resolution of the temperature	0.1
	6.5 5.7		0.1
Usable immersion depth in. Weight lbs		Setting the resolution of the temperature	0.1 +20 +150
Usable immersion depth in.	5.7	Setting the resolution of the temperature display °C	
Usable immersion depth in. Weight lbs	5.7	Setting the resolution of the temperature display °C Working temperature range °C	+20 +150

For applications near or below ambient temperature: use a cooling coil or JULABO immersion cooler.

All Benefits



100% Checked.

100% testing. 100% quality. Each JULABO Circulator undergoes thorough quality testing before leaving the factory.



Green technology. Development consistently applied environmentally friendly materials and technologies.



JULABO. Quality. Highest standards of quality for a long product life.



Satisfied customers. 11 subsidiaries and more than 100 partners worldwide guarantee fast and qualified JULABO support.



Quick start. Individual JULABO consultation and comprehensive manuals at your disposal.



ATC. Absolute Temperature Calibration, 1-point calibration (CD).



below the display. Easily change between internal and external circulation.

Internal. External.



Modern. Reliable.

High-grade components in every JULABO Circulator - platinum sensors, proven motor technology, CAN-Bus communication and much more.

The pump is controlled via a lever located directly

Services 24/7.

Around the clock availability. You can find suitable accessories, data sheets, manuals, case studies, and more at www.julabo.us.



Quick support

If an error occurs, the integrated Black-Box function permits fast diagnosis by the JULABO service team



2





Timer. Integrated.

CORIO circulators include an integrated timer function. When the set time has elapsed, a signal sounds and the device switches off. Setting range: 0 ... 999 minutes.



Connection. Easy.

Inclined pump connections (M16×1) facilitate the connection of applications. Each unit includes 2 barbed fittings of 8/12 mm diameter each.



Exact.

You can rely on it. PID1 control and 'Active Cooling Control' make the new CORIO precise and perfect.



Switch on. And off you go.

Intelligent operating concept. Ready for operation with just a few quick and easy steps.



Connectivity.

Remote control made easy. CORIO CD circulators feature a USB connection.



Brilliant.

Very bright display makes it easy to read even from a distance.



Everything at the front.

All operating controls and safety functions are accessed easily and comfortably from the front.



Locked in. The lockable power plug guarantees safe connection. More process safety.

Early warr Maximum

Early warning system for low liquid level Maximum safety for applications, optical and audible alarm, allows user to refill bath fluid before the unit shuts down