Julabo

FW95-SL Ultra-Low Refrigerated / Heating Circulator

Ultra-Low Refrigerated Circulators with SL circulator impress with their high heating, cooling and pump capacities for external temperature control applications.

Product features

- VFD COMFORT DISPLAY
- Keypad for setpoints, warning/safety values and menu functions
- ICC (Intelligent Cascade Control), self-optimizing temperature control
- TCF Temperature Control Features to optimize the control behavior
- ATC3 3-Point-Calibration
- Pt100 External sensor connection for measurement and control
- SMART PUMP, electronically adjustable pump stages
- Adjustable high temperature cut-out, visible via display
- RS232/RS485 interface for online communication
- Integrated programmer for 6 x 60 program steps
- Connections for solenoid valve and HSP booster pump
- Proportional cooling control



Performance values

208V/3PPE/60Hz (Without Plug)	
Heating capacity kW	3
Pump capacity flow rate l/min	22 26
Pump capacity flow pressure psi	5.8 10.2
Maximum suction psi	2.9 5.8
Current Draw A	28

Order No.		9352796.16.N			
Cooling capacity (Ethanol)					
°C	0	-20	-40	-60	-80
kW	1.7	1.5	1.3	1.1	0.36

*Performance specifications measured in accordance with DIN12876. Cooling capacities up to 20 °C measured with ethanol; over 20 °C with thermal oil unless otherwise specified. Performance specifications apply at an ambient temperature of 20 °C. Performance values may differ with other bath fluids.

Refrigerant stage 1		Refrigerant stage 2	
Refrigerant	R404A	Refrigerant	R23
Filling volume g	1000	Filling volume g	380
Global Warming Potential for R404A	3922	Global Warming Potential for R23	14800
Carbon dioxide equivalent	t 3.922	Carbon dioxide equivalent	t 5.624



Technical data

Available voltage vers	ions	Cooling	
Order No.	9 352 796.N	Cooling of compressor	2-stage Water
Available voltage versions:		Cooling water temperature max. °C	20
9352796.07.N	400V/3PNPE/50Hz (Plug 32A CEE)	Cooling water pressure max. psi	87
(R404A)	Cooling water difference pressure psi	50.8	
9352796.16.N	208V/3PPE/60Hz (Without Plug) (R404A)		

Bath		Other	
Bath cover	integrated	Classification	Classification III (FL)
		IP Code	IP 21
Electronics		Dimensions and volumes	
Digital interface	Profibus optional	Weight lbs	443.1
External pt100 sensor connection	integrated	Dimensions in. ($W \times L \times H$)	23.2 x 29.9 x 45.7
Integrated programmer	6x60 steps	Filling volume l	22
Temperature control	ICC	Pump connections	M16x1 male
Absolute temperature calibration	3 Point Calibration		
Temperature display	VFD		
Temperature setting	Keypad		
Temperature values		Included in delivery	
Setting the resolution of the temperature display °C	0.01	2 Barbed fittings for tubing 8 and 12 mm II M16x1 male)	D. (Pump connections
Working temperature range °C	-95 0		
Temperature stability °C	±0.05		
Ambient temperature °C	+5 +40		
Temperature display resolution °C	0.01		

All Benefits



100% Checked.

100% testing. 100% quality. Each JULABO Circulator undergoes thorough quality testing before leaving the factory.



JULABO. Quality. Highest standards of quality for a long product life.



Quick start.

Development consistently applied environmentally

Individual JULABO consultation and comprehensive manuals at your disposal.

friendly materials and technologies.

Green technology.





Satisfied customers.

11 subsidiaries and more than 100 partners worldwide guarantee fast and qualified JULABO support.



Intelligent temperature control. Intelligent cascade control - automatic and self-

optimizing adaptation of the PID control parameters with external stability of +/- 0.05 °C.



Early warning system for high/low temperature limits

Maximum safety for applications, optical and audible alarm, convertible to automated cut-off function



For flammable bath fluid

Classification III (FL) according to DIN 12876-1



100 % Cooling capacity

'Active Cooling Control' for cooling available throughout the entire working temperature range, fast cool-down even at higher temperatures



Condensation and ice protection A heated cover plate prevents condensation

A heated cover plate prevents condensation or ice build-up in the bath



ATC3. Calibration.

'Absolute Temperature Calibration' for compensating a physically caused temperature difference, 3-point calibration.



Services 24/7. Around the clock availability. You can find suitable accessories, data sheets, manuals, case studies, and more at www.julabo.us.



Connection of additional equipment Stakei connections for solenoid valve, HSP booster pump and HST booster heater



Clever pump system Reliable and consistent pump capacity, electronically adjustable pump stages



Process. Under control. Full control of the dynamic, access to all important control parameters for individual process optimization.



Energy saving cooling

Proportional cooling control for automatic adjustment of cooling power or temporary switchoff of compressor as needed to save up to 90 % energy in comparison to unregulated cooling machines



Control of the external application

External Pt100 sensor connection for precise measurement and control directly in the external application

