# Julabo

## MAGIO MS-1000FW Refrigerated / heating circulator

As with all circulators from the MAGIO range, the refrigerated circulators stand out thanks to their premium quality, high performance and intuitive operation. The devices offer extra strong pressure and suction pumps, thus fulfilling the highest demands for temperature control of external applications. Whether in basic research, material testing or technical systems – the MAGIO refrigerated circulators offer high-tech solutions for high customer requirements.

#### High resolution TFT touch display

The modern TFT touch display gives you all the important information at a glance. Three large, predefined main screens clearly display data and graphics with various application priorities. Menu navigation is self-explanatory, arranged by relevance to daily operations and easy to operate with the touch of a finger. The in-built help function provides detailed support in case of additional questions.





#### **Product features**

- Ideal for demanding external applications
- Simple control of complex applications
- Continuously adjustable, extremely powerful pressure / suction pump
- Flow rate 16 ... 31 I / min, pressure 0.24 ... 0.92 bar, suction 0.03 ... 0.4 bar
- Large, high-resolution TFT touch display with multilingual user interface
- Parts being in contact with the medium made of stainless steel
- Integrated programmer
- Integrated external Pt100 connection
- USB connection
- RS232 interface for online communication
- Ethernet
- analog interfaces (accessory)
- Class III (FL) according to DIN 12876-1
- Modbus
- Profibus DP (Accessory)
- RS232/RS485 interface for online communication
- Connections for solenoid valve



#### Performance values

115V/60Hz (Nema N5-20 Plug)	
Heating capacity kW	1
Viscosity max. cST	70
Pump capacity flow rate I/min	16 31
Pump capacity flow pressure psi	3.5 13.3
Maximum suction psi	-0.45.8
Power A	16

Order No.				9032727	7.02		
Cooling cap	oacity (Etha	nol)					
°C	20	10	0	-10	-20	-30	-40
kW	1	0.96	0.96	0.7	0.51	0.25	0.11

\*Performance specifications measured in accordance with DIN12876. Cooling capacities up to 20 °C measured with ethanol; over 20 °C with thermal oil unless otherwise specified. Performance specifications apply at an ambient temperature of 20 °C. Performance values may differ with other bath fluids.

Refrigerant stage 1	
Refrigerant	R449A
Filling volume g	190
Global Warming Potential for R449A	1397
Carbon dioxide equivalent t	0.265

#### **Technical data**

Available voltage ve	ersions	Bath	
Order No.	9 032 727	Bath tank	Stainless steel
Available voltage version	ns:	Bath cover	integrated
9032727.02	115V/60Hz (Nema N5-20 Plug) (R449A)	Usable bath opening in. (W x L / D)	7.1 x 5.1 / 5.9
9032727.05	200-230V/50-60Hz (CH Plug Type SEV 1011) (R449A)		
9032727.04	200-230V/50-60Hz (UK Plug Type BS1363A) (R449A)		
9032727.33	200-230V/50-60Hz (Schuko Plug - CEE 7/4 Plug Type F) (R449A)		
9032727.33.chn	200-230V/50-60Hz (CN Plug) (R449A)		

Cooling		Other	
Cooling of compressor	1-stage Water	Classification	Classification III (FL)
Cooling water temperature max. °C	30	IP Code	IP 21
Cooling water pressure max. psi	87	Pump function	Pressure Suction Pump
Cooling water difference pressure psi	4.4	Pump type	Immersion Pump
1.0Cooling water consumption l/min	1.1		

## Julabo

User Interface Language

Chinese, Czech, Dutch, English, French, German, Italian, Japanese, Korean, Polish, Portuguese, Russian, Spanish, Turkish

Electronics	
Interfaces	Ethernet, Modbus, RS232, RS485, USB
External pt100 sensor connection	integrated
Integrated programmer	8x60 steps
Temperature control	ICC
Absolute temperature calibration	10 Point Calibration
Temperature displayTemperature display	7" TFT Touchscreen
Temperature settingTemperature setting	Touchscreen
Electronic Timer hr:min	00:00 99:59

Weight Ibs	119.3
Dimensions in. ( $W \times L \times H$ )	16.5 x 19.3 x 29.1
Filling volume I	5 7.5
Pump connections	M16x1 male

Temperature values	
Setting the resolution of the temperature display $^\circ\mathrm{C}$	0.01
Working temperature range °C	-50 +200
Temperature stability °C	±0.01
Ambient temperature °C	+10 +40
Temperature display resolution °C	0.01

#### Included in delivery

2 Barbed fittings for tubing 8 and 12 mm ID. (Pump connections M16x1 male)

#### All Benefits



#### Quick support

If an error occurs, the integrated Black-Box function permits fast diagnosis by the JULABO service team



### Many interfaces.

Straight-forward remote control, data management, and integration into process structures. USB, Ethernet, RS232, SD card, and alarm off are permanently integrated. Further interfaces available as accessories.



#### Maximum safety.

Classification III according to DIN12876-1 enables safe operation, even with flammable fluids. Automatic switch-off in the event of high temperature or low liquid level.



Multi-lingual. Operation in multiple languages.



#### Intelligent temperature control.

Intelligent cascade control - automatic and selfoptimizing adaptation of the PID control parameters with external stability of +/- 0.05 °C.



#### Touch display. Perfect operation.

With the touch display, the user always has an overview of all values and functions. The intuitive and multilingual menu structure enables perfect control.

#### Space saving. Free up space.



Place your JULABO Circulator right next to an application, another unit, or wall. That saves space. This is made possible by eliminating vents and connections on the sides.

#### Programmer. Integrated.

The integrated internal programmer makes it possible to automatically run temperature time profiles.

application.





**Temperature. Under control.** External Pt100 sensor connection for precise



Analog I/O. Analog interfaces for integration into process control systems (optional).

Early warning - visual and acoustic - of critical

states increases process stability.

Everything made of stainless steel.

entirely made of stainless steel.



Fill level. Monitored. Fill level indicator on the display for heat-transfer liquid.

measurement and control directly in the external



**Process. Under control.** Full control of the dynamic, access to all important control parameters for individual process optimization.



Energy-saving. The high-quality insulation of all relevant components saves energy.



#### Wide range.

Refrigerated and heating circulator in various combinations, circulator in various sizes. Maximum flexibility through a large selection of accessories.



#### Most powerful pump.

The integrated pressure/suction pump with performance values of 0.9 bar and -0.4 bar is the most powerful in its class and continuously adjustable.



#### Highest measuring accuracy

'Absolute Temperature Calibration' for manual compensation of a temperature difference, 10point calibration



#### Connection. Easy.

Process stability.

Stable. Mobile.

Inclined pump connections (M16×1) facilitate the connection of applications. Each unit includes 2 barbed fittings of 8/12 mm diameter each.

Quality and material compatibility at the highest

level. All parts in contact with the medium are



#### Condensation protection.

Superb design solution. Integrated ventilation directs air over the bath lid and minimizes condensation.

