MAGIO MX-BC12 HST Heating Circulator

MAGIO heating circulators with booster heater are characterized by premium quality and high performance combined with intuitive usability. The MAGIO offers powerful pressure suction pumps to meet the most challenging demands for temperature control of external applications.

High Resolution TFT Touch Display

The color TFT touch display gives you all the important information immediately. Three large, predefined main screens clearly display data and graphics with various application priorities. Menu navigation is self-explanatory, arranged by relevance to daily operations and easy to operate with the touch of a finger. The built-in help function provides detailed support in case of additional questions.



Product Features

- For external applications only
- Continuously adjustable, powerful pressure / suction pump
- · Large, high-resolution TFT touch display with multilingual user interface
- All wetted components made from stainless-steel
- Flow rate 16 ... 31 L / min | pressure 3.5 ... 13.3 psi, suction -2.3 ... -5.8 psi
- Integrated programmer
- Integrated external Pt100 connection
- USB B connection for PC integration | USB A connection for data-logging
- RS232 / RS485 interface for online communication
- Ethernet, Modbus TCP/IP
- Optional: Analog interfaces (accessory p/n 9900100)
- High-quality thermal insulation of the bath tanks
- Very quiet
- Easy operation
- Strong heating capacity of 9kW at 230V or 7.4kW at 208V
- Connections for solenoid valve
- Integrated pump connection M16×1

Performance values

| 200-230V/3PPE/50-60Hz (60Hz shipped bare wire, MX plugs into HST control box) | | | |
|---|-----------------------------|--|--|
| Heating capacity kW | 9kW at 230V 7.4kW at 208V | | |
| Viscosity max. cSt | 70 | | |
| Pump capacity flow rate I/min | 16 31 | | |
| Pump capacity flow pressure psi | 3.5 13.3 | | |
| Maximum suction psi | -0.45.8 | | |
| Power consumption A | 28A | | |

Technical data

| Available voltage vers | ions | Bath | |
|---|-------------------------------------|--|---|
| Order No. 9033512.16.HST 9033512.16.HST.CSA/UL (with CSA C22.2 and UL 61010-1 NRTL approvals) | Bath tank | Stainless steel | |
| | Usable bath opening in. (W x L / D) | 8.7 x 5.9 / 7.9 (without HST installed) | |
| Available voltage versions: | | | , |
| 9033512.16.HST | 200-230V/3PPE/60Hz (Bare wire) | | |

| Cooling | | Other | |
|-----------------------|------------|-------------------------|--|
| Cooling coil integrat | integrated | Classification | Classification III (FL) |
| | | IP Code | IP 21 |
| | | Pump function | Pressure/Suction Pump |
| | | Pump type | Immersion Pump |
| | | User Interface Language | Chinese, Czech, Dutch, English, French, German, Italian, Japanese, Korean, Polish, Portuguese, Russian, Spanish, Turkish |

| Electronics | | Dimensions and volumes | |
|----------------------------------|--|--|-------------------------------------|
| Interfaces | erfaces Ethernet, Modbus TCP/IP, RS232, RS485, Stakei, USB | Weight lbs | 50 |
| | | Dimensions in. (W × L × H) MX-BC12 w/HST in bath opening HST control box | 13 x 19.3 x 18.5 7.9 x 6.7 x 7.5 |
| | | Filling volume l | 8.5 12 |
| External pt100 sensor connection | integrated | | |
| Integrated programmer | 8x60 steps | - | |
| Temperature control | ICC | Pump connections | M16x1 male |
| Absolute temperature calibration | 10 Point Calibration | | |
| Temperature display | 7" TFT Touchscreen | | |
| Temperature setting | Touchscreen | | |
| | 00:00 99:59 | | |



| Temperature values | |
|--|----------|
| Setting the resolution of the temperature display °C | 0.01 |
| Working temperature range °C | +20 +300 |
| Temperature stability °C | ±0.01 |
| Ambient temperature °C | +10 +40 |
| Temperature display resolution °C | 0.01 |

Included in delivery

2 Barbed fittings for tubing 8 and 12 mm ID. (Pump connections M16x1 male)

All Benefits



Intelligent temperature control. Intelligent cascade control - automatic and selfoptimizing adaptation of the PID control parameters with external stability of +/- 0.05 °C.



Many interfaces.

Straight-forward remote control, data management, and integration into process structures. USB, Ethernet, RS232, SD card, and alarm off are permanently integrated. Further interfaces available as accessories.

Touch display. Perfect operation.

With the touch display, the user always has an overview of all values and functions. The intuitive and multilingual menu structure enables perfect control.



Space saving. Free up space.

Place your JULABO Circulator right next to an application, another unit, or wall. That saves space. This is made possible by eliminating vents and connections on the sides.



Programmer. Integrated.

The integrated internal programmer makes it possible to automatically run temperature time profiles.



Temperature. Under control. External Pt100 sensor connection for precise measurement and control directly in the external application.



Process stability.

Early warning - visual and acoustic - of critical states increases process stability.



Maximum safety.

Classification III according to DIN12876-1 enables safe operation, even with flammable fluids. Automatic switch-off in the event of high temperature or low liquid level.



Multi-lingual.

Operation in multiple languages.



Analog I/O. Analog interfaces for integration into process control systems (optional).



Fill level. Monitored. Fill level indicator on the display for heat-transfer liquid.



Process. Under control. Full control of the dynamic, access to all important control parameters for individual process optimization.





Stable. Mobile.



Everything made of stainless steel. Quality and material compatibility at the highest level. All parts in contact with the medium are entirely made of stainless steel.



Wide range.

Refrigerated and heating circulator in various combinations, circulator in various sizes. Maximum flexibility through a large selection of accessories.



Connection. Easy.

Inclined pump connections (M16×1) facilitate the connection of applications. Each unit includes 2 barbed fittings of 8/12 mm diameter each.



100% Checked.

100% testing. 100% quality. Each JULABO Circulator undergoes thorough quality testing before leaving the factory.



JULABO. Quality.

Highest standards of quality for a long product life.



Satisfied customers.

11 subsidiaries and more than 100 partners worldwide guarantee fast and qualified JULABO support.



Energy-saving.

The high-quality insulation of all relevant components saves energy.



Most powerful pump.

The integrated pressure/suction pump with performance values of 0.9 bar and -0.4 bar is the most powerful in its class and continuously adjustable.



Condensation protection.

Superb design solution. Integrated ventilation directs air over the bath lid and minimizes condensation.



Highest measuring accuracy 'Absolute Temperature Calibration' for manual compensation of a temperature difference, 10point calibration



Green technology.

Development consistently applied environmentally friendly materials and technologies.



Ouick start.

Individual JULABO consultation and comprehensive manuals at your disposal.



Services 24/7.

Around the clock availability. You can find suitable accessories, data sheets, manuals, case studies, and more at www.julabo.com.