

CORIO™

CD PRODUCT LINE



CORIO CD

SKU 9012000



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **to 50**
Dimensions **5.2 x 6.5 x 14.2**

CD-1000F

SKU 9012707



Working Temp Range (°C) **-40... +150**
Cooling Capacity **1 kW @ 20°C**
1 kW @ 10°C | 0.98 kW @ 0°C
0.75 kW @ -10°C | 0.53 kW @ -20°C
0.27 kW @ -30°C | 0.13 kW @ -40°C
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **5... 7.5**
Dimensions **16.5 x 19.3 x 27.6**

CD-1001F

SKU 9012708



Working Temp Range (°C) **-38... +100**
Cooling Capacity **1 kW @ 20°C**
0.95 kW @ 10°C | 0.9 kW @ 0°C
0.63 kW @ -10°C | 0.35 kW @ -20°C
0.13 kW @ -30°C
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **42... 56**
Dimensions **17.7 x 25.2 x 37.4**

CD-200F

SKU 9012701



Working Temp Range (°C) **-20... +150**
Cooling Capacity **0.22 kW @ 20°C**
0.17 kW @ 0°C | 0.13 kW @ -10°C
0.06 kW @ -20°C
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **3... 4**
Dimensions **9.1 x 15.4 x 25.6**

CD-201F

SKU 9012702



Working Temp Range (°C) **-20... +150**
Cooling Capacity **0.22 kW @ 20°C**
0.16 kW @ 0°C | 0.13 kW @ -10°C
0.06 kW @ -20°C
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **3... 4**
Dimensions **17.4 x 16.2 x 17.4**

CD-300F

SKU 9012703



Working Temp Range (°C) **-25... +150**
Cooling Capacity **0.31 kW @ 20°C**
0.28 kW @ 0°C | 0.2 kW @ -10 °C
0.11 kW @ -20°C
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **3... 4**
Dimensions **9.5 x 16.6 x 26**

CD-600F

SKU 9012704



Working Temp Range (°C) **-35... +150**
Cooling Capacity **0.6 kW @ 20°C**
0.46 kW @ 0°C | 0.26 kW @ -10 °C
0.18 kW @ -20°C | 0.06 kW @ -30 °C
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **5... 7.5**
Dimensions **13 x 18.5 x 27.2**

CD-601F

SKU 9012705



Working Temp Range (°C) **-40... +150**
Cooling Capacity **0.6 kW @ 20°C**
0.55 kW @ 10 °C | 0.4 kW @ 0°C
0.29 kW @ -10°C | 0.18 kW @ -20°C
0.06 kW @ -30°C | 0.01 kW @ -40°C
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **8... 10**
Dimensions **14.2 x 18.1 x 29.1**

CD-900F

SKU 9012706



Working Temp Range (°C) **-40... +150**
Cooling Capacity **0.9 kW @ 20°C**
0.9 kW @ 10 °C | 0.8 kW @ 0°C
0.55 kW @ -10°C | 0.35 kW @ -20°C
0.15 kW @ -30°C | 0.02 kW @ -40°C
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **21... 30**
Dimensions **15.4 x 24.4 x 29.5**

CD-B13

SKU 9012413



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **9... 13**
Dimensions **15 x 15.8 x 16.6**

CD-B17

SKU 9012417



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **13... 17**
Dimensions **15 x 15.8 x 18.6**

CD-B19

SKU 9012419



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **14... 19**
Dimensions **15 x 23 x 16.6**

CORIO™

CD PRODUCT LINE



CD-B27

SKU 9012427



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **17... 27**
Dimensions **15 x 23 x 18.6**

CD-B33

SKU 9012433



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **26... 39**
Dimensions **36 x 14.2 x 17**

CD-B39

SKU 9012439



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **35... 41**
Dimensions **21.3 x 13.4 x 22.5**

CD-B5

SKU 9012405



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **3.5... 5**
Dimensions **9.1 x 15 x 16.2**

CD-BC12

SKU 9012512



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **8.5... 12**
Dimensions **13 x 19.3 x 18.5**

CD-BC26

SKU 9012526



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **19... 26**
Dimensions **15.4 x 24.5 x 18.9**

CD-BC4

SKU 9012504



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **3... 4.5**
Dimensions **9.1 x 16.2 x 16.6**

CD-BC6


SKU 9012506



Working Temp Range (°C) **+20... +150**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **4.5... 6**
Dimensions **9.5 x 17.4 x 18.5**

CD-BT19

SKU 9012319



Working Temp Range (°C) **+20... +100**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **14... 19**
Dimensions **15 x 22.9 x 15**

CD-BT27

SKU 9012327



Working Temp Range (°C) **+20... +100**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **20... 27**
Dimensions **15 x 22.9 x 17**

CD-BT5

SKU 9012305



Working Temp Range (°C) **+20... +100**
Cooling Capacity **N/A**
Heating Capacity (kW) **1**
Pump Capacity Flow Rate (l/min) **15**
Pump Capacity Flow Pressure (psi) **5.1**
Temperature Stability (°C) **±0.03**
Filling Volume Liters **3.5... 5**
Dimensions **9.1 x 15 x 15**